




MetaSinglePort & WR


Disposable Laparoscopic Single Port Access Device Filter
& Wound Retractor



Product Overview

MetaSinglePort is a product used for a single-hole laparoscopic surgery. Its specially designed locking plate mechanism is simple to use and it maximizes the user's vision of surgery by using clear & transparent window. Advanced filter system has been developed for moisture adsorption and deodorization.

MetaSinglePort enables operations With minimal invasion surgery



Key Features

- Easy to use locking plate
- Clear & Transparent window
- Advanced filter (large capacity, long-term usability)



Filter Performance Test

High ammonia deodorization rate and TVOC adsorption

Product name	MetaSinglePort	
Test	Ammonia	TVOC 총휘발성유기화합물
Time (hours)	Removal efficiency	Deodorization rate
3 hours	98%	97%
4 hours	99%	98%
5 hours	99.8%	99%

* Table 1. Test results based on company standards.

Order Information

Model No	Description	Specification	Quantity
MSP01A	4Channel Port	5-10mm 3channel	5EA / 1BOX
		10-12mm 1channel	
	Wound Retractor	ID: 60mm, L: 180mm	
	Filter	-	
MWR01A	Wound Retractor	ID: 60mm, L: 180mm	10EA / 1BOX

meta Surgical

Head Office

Meta Biomed Co.,Ltd.

270 Osongsaeangmyeong1-ro,
Osong-eup, Heungdeok-gu,
Cheongju-si, Chungbuk, Korea
Phone: +82 43 230 8846
Fax: + 82 43 230 8969
boram@meta-biomed.com
www.meta-biomed.com

Meta Biomed Inc.

121 Independence Lane
Chalfont, PA 18914, USA
Phone: +1 267 282 5893
america@metabiomed-inc.com
www.buymetabiomedonline.com

Meta Biomed Europe

Haus der Wirtschaft, Raum 3.11
Wiesenstraße 35
45473 Mülheim an der Ruhr, Germany
Phone: +49 208 309 9190
europe@metabiomed-inc.com
www.buymetabiomedonline.com